BTSID1 Product Change Notifications

Product change notifications (PCN) is the way of informing customers of upcoming changes related to Teltonika products. This includes outward appearance adjustments, hardware improvements, firmware changes, changes to package contents, etc. This page stores information regarding changes to the BTSID1 tracking device. Below you will find a list of both upcoming and already implemented changes sorted from newest to oldest (top to bottom).

Contents

- 1 2024.02.22: Base firmware release
- 2 2023.11.20: Hall sensor added to base BTSID1 hardware
- 3 2023.11.16: Base firmware release
- 4 2023.05.01: Base firmware release
- 5 2022.12.14: Base firmware update
- 6 2022.09.09: BOT Laser Printing Changes
- 7 2022.07.21: Base firmware update
- 8 2022.07.21: Advanced configuration possibilities in manufacturing
- 9 2022.01.28: BOT Laser Printing Changes

2024.02.22: Base firmware release

We would like to inform you about the base firmware release for BTSX devices starting from February 23rd, 2024

New production firmware version will be placed in manufacturing, version 1.2.9.R.9 will be used for new order codes with hibernate feature.

Change description

Change description		
Change type	Firmware update	
Detailed description	Changes from last production firmware are listed in <u>BTS firmware</u> <u>errata</u> page	
	Old version:	New version:
	BTSX.Ver.1.2.10.R.04	1.2.9.R.9
Change reasons	BTSX.Ver.1.2.10.R.04 found to have issue with hibernate	
Product release date	2023 02 23	
Order codes with new firmware	BTSMP1SAMPLE1	SPQ:1 sleep (envelope)
	BTSMP1SAMPLE2	SPQ:1 + magnet (box)
	BTSMP14NB802	SPQ:1 + magnet (box)
	BTSMP14NE502	SPQ:10
	BTSMP14NM402	SPQ:100

2023.11.20: Hall sensor added to base BTSID1 hardware

We would like to inform you about the base hardware update for BTSID1 devices starting from November 30th, 2023

Updated hardware version will enable use of device magnetic field sensing functions. Like hibernate or door status detection.

Change description

Change type Hardware update

Old version: BTSID1 with no New version: BTSID1 with hall Detailed description

sensors

1. Allow all devices to have Magnetic sensing capabilities for

Change reasons hibernate function

2. Simplify manufacturing by having only one BTSID1 version.

Change date 2023 11 30

Order codes with new

All codes without hibernate will have analogs created. Check for analog code with your Sales manager hardware

Change impact

No risk factors encountered with this change Risk assessment

Suggested implementation plan No implementation plan needed

Acknowledgement of PCN receipt

If no feedback is received within two weeks after the issue date of this notification - Teltonika may accept that this change has been tacitly accepted and can implement the change as indicated above

2023.11.16: Base firmware release

We would like to inform you about the base firmware release for BTSX devices starting from November 30th, 2023

New production firmware version will be placed in manufacturing, version 1.2.10.R.4 will be used for new order codes with hibernate feature.

Change description

Change type Firmware update

Improvements from last production firmware are listed in BTS firmware

errata page

New version: Old version:1.2.9.R.8

BTSX.Ver.1.2.10.R.04

New base version in factory Change reasons

Product release date 2023 11 30

Detailed description

Order codes with new

BTSID1SAMPLE1 SPQ:1 sleep (envelope) BTSID1SAMPLE2 SPQ:1 + magnet (box)BTSID17RB802 SPQ:1 + magnet (box)

firmware BTSID17RE502 SPQ:10 BTSID17RM402 SPO:100

2023.05.01: Base firmware release

We would like to inform you about the base firmware release for BTSX devices starting from May 1st, 2023.

New production firmware version will be placed in manufacturing, version 1.2.9.R.8 will be used for new order codes to support hibernate feature.

Change description

Change type Firmware update

Improvements from last production firmware are listed in **BTS**

<u>firmware errata</u> page

Detailed description

Old version:

BTSX.Ver.1.2.9.R.8

Change reasons New base version in factory

Product release date 2023 May

BTSID1SAMPLE1 SPQ:1 sleep (envelope)
BTSID1SAMPLE2 SPQ:1 + magnet (box)
BTSID17RB802 SPQ:1 + magnet (box)

Order codes with new firmware

BTSID17RE502 SPQ:10 BTSID17RM402 SPQ:100

2022.12.14: Base firmware update

We would like to inform you about the base firmware update for BTSX devices starting from January 1st, 2023. New production firmware version will be 1.2.8.R.2

Change description

Change type Firmware update

Improvements from last production firmware are listed in BTS firmware

x x

errata page

Detailed description

Old version: New version: BTSX.Ver.1.2.8.R.1 BTSX.Ver.1.2.8.R.2

Change reasons New base version in factory

Product release date 2023 January

2022.09.09: BOT Laser Printing Changes

We would like to inform you, that laser printing of the BTSID1 device was updated in the factory in September 2022.

Change description

Change type External/visual

Bot laser printing design changed

Detailed description

Old version: New version:

x x

1. Added SN number to casing TOP

Change reasons

1. Added Tales continued to casing for

2. Added Telec certification to casing BOT

Product release

date October 2022

Affected products

Model Affected order codes Order code changes
BTSID1 BTSID1*****

New laser engraving

Change impact

Risk assessment No risk factors encountered with this change

Suggested Due to changes in device appearance we strongly recommend to contact your customs agent with our change log information, to avoid any unexpected delays

plan in customs clearance process

Acknowledgement of PCN receipt

If no feedback is received within two weeks after the issue date of this notification - Teltonika may accept that this change has been tacitly accepted and can implement the change as indicated above

2022.07.21: Base firmware update

We would like to inform you about the base firmware update for BTSX devices starting from July 21st, 2022. New production firmware version will be 1.2.8

Change description

Change type Firmware update

Key improvements from last production firmware are listed in BTS firmware

errata page

Detailed description Old version: New version:

BTSX.Ver.1.2.1 BTSX.Ver.1.2.8

Change reasons New base version

Product release date 2022 July

Affected products

Change type Default order code

BTSID1 BTSID1******* order codes with base Firmware

Old version: New version:

Detailed description BTSID14KB801 BTSID1W1B801

BTSID14KE501 BTSID1W1E501

Change reasons
Important default configuration parameters (advertising interval) changed

2022.07.21: Advanced configuration possibilities in manufacturing

We would like to inform you that starting from July 21st, 2022. Manufacturing is capable of advanced pre-configuration options for BTSXXX devices.

It's now possible to configurate individual parameter values according to clients' demands, e.g.:

1. It's possible to write MAC into any other information advertising parameter that fits 6 bytes of information for example you can include MAC into UUID of the beacon into of name of the beacon.

2. Also during production values can be generated from specific value and according to your logic. For example, you can generate minor from 0000 to FFFF.

2022.01.28: BOT Laser Printing Changes

We would like to inform you, that laser printing of the BTSID1 device was updated in the factory in January 2022.

Change description

Change type External/visual

Bot laser printing design changed

Detailed description Old version:

New version:

1. Added Anatel certificate logo and number.

Change reasons 2. Added E-MARK certificate logo and number.

3. Added FCC ID.

Product release

date January 2022

Affected products

Model Affected order codes Order code changes
BTSID1 BTSID1*****

New laser engraving

Change impact

Risk assessment No risk factors encountered with this change

Suggested Due to changes in device appearance we strongly recommend to contact your implementation customs agent with our change log information, to avoid any unexpected delays

plan in customs clearance process

Acknowledgement of PCN receipt

If no feedback is received within two weeks after the issue date of this notification - Teltonika may accept that this change has been tacitly accepted and can implement the change as indicated above